

tri▼temp  
new!



### P-2420

#### OVERALL SIZE

30"W x 21¼"D x 25¼"H

#### COMPARTMENT SIZE

- 4.75 cubic ft.
- 26½"W x 17"D x 18¼"H

#### • tri▼temp Range:

- Irrigation: 90°F - 150°F
- Injection: 90°F - 104°F
- Blanket: 90°F - 160°F

#### • Total Capacity:

- 20 1-liter bags or
- 24 1-liter bottles
- 9 - 12 blankets

#### • Data Logger

#### • Caster (Option)

tri▼temp  
new!



### P-2430

#### OVERALL SIZE

30"W x 27½"D x 25¼"H

#### COMPARTMENT SIZE

- 5.5 cubic ft.
- 26½"W x 19¾"D x 18¼"H

#### • tri▼temp Range:

- Irrigation: 90°F - 150°F
- Injection: 90°F - 104°F
- Blanket: 90°F - 160°F

#### • Total Capacity:

- 30 1-liter bags or
- 48 1-liter bottles
- 10 - 14 blankets

#### • Data Logger

#### • Caster (Option)

## warming cabinets

blanket & fluid

#### features

Combination warmers feature tri-temp chamber allowing for 3 different modes of warming, LED display, window doors, and control lock programming.



### P-2010-S

#### OVERALL SIZE

16½"W x 25¼"D x 21½"H

#### COMPARTMENT SIZE

- 2.3 cubic ft.
- 13¼"W x 21¾"D x 14"H

#### • Temp Range:

- Blankets: 98°F - 200°F

#### • Total Capacity:

- 5 - 7 blankets

#### • Solid door

#### • Caster (Option)



### P-2012

#### OVERALL SIZE

18½"W x 22"D x 22"H

#### COMPARTMENT SIZE

- 2.5 cubic ft.
- 15½"W x 18½"D x 15"H

#### • Temp Range:

- Blanket: 90°F - 160°F

#### • Total Capacity:

- 4 - 6 blankets

#### • Caster (Option)

Heated center  
shelf



### P-2022

#### OVERALL SIZE

18½"W x 22"D x 28"H

18½"W x 23"D x 34½"H (with casters)

#### COMPARTMENT SIZE

- 3.5 cubic ft.
- 15½"W x 18½"D x 21"H\*

\* (9 9/16" H above and below the center shelf)

#### • Temp Range:

- Blanket: 90°F - 160°F

#### • Total Capacity:

- 10 - 12 blankets

#### • Caster (Option)

shown with left door  
hidge mount option



### P-2032

#### OVERALL SIZE

24"W x 24½"D x 35½"H

24"W x 25¼"D x 42"H (with casters)

#### COMPARTMENT SIZE

- 7.5 cubic ft.
- 21"W x 21"D x 28"H

\* (13 3/16" H above and below the center shelf)

#### • Temp Range:

- Blanket: 90°F - 160°F

#### • Total Capacity:

- 20 - 25 blankets

#### • Caster (Option)



## P-2040

### OVERALL SIZE *(without bumper)*

24"W x 34 $\frac{1}{4}$ "D x 72 $\frac{3}{4}$ "H

### COMPARTMENT SIZE

- 7.7 cubic ft. upper
- 19"W x 26 $\frac{1}{2}$ "D x 26 $\frac{3}{4}$ "H
- 7.7 cubic ft. lower
- 19"W x 26 $\frac{1}{4}$ "D x 26 $\frac{3}{4}$ "H

- **Temp Range:**
  - Blankets: 98°F - 160°F
- **Total Capacity:**
  - 40-50 blankets
- Two compartments, each with individual controls
- Auto time programming
- 5" casters
- Pass Through Option Available

# warming cabinets

blanket & fluid

## features

Combination fluid and blanket warmers feature LED display, window doors, and control lock programming.

- No fans or lint traps
- Two compartments, each with individual controls



## P-2147

### OVERALL SIZE

- 30"W x 28 $\frac{1}{2}$ "D x 75"H (caster standard)
- 30"W x 28 $\frac{1}{2}$ "D x 68"H (feet - option)

### COMPARTMENT SIZE

- 5.5 cu. ft. fluid compartment
- 26 $\frac{1}{2}$ "W x 23 $\frac{3}{4}$ "D x 15"H
- 13.25 cu. ft. blanket compartment
- 26 $\frac{3}{4}$ "W x 23 $\frac{3}{4}$ "D x 36"H

- **tri▼temp Range:**
  - Blankets: 90°F - 160°F
  - Irrigation: 90°F - 150°F
  - Injection: 90°F - 104°F
- **Total Capacity:**
  - 24-36 1-liter bags
  - 30 1-liter bottles
  - 10-14 blankets (upper chamber)
  - 30-40 blankets (lower chamber)
- Separate **tri▼temp** temperature controls for each compartment
- 5" casters
- Data Logger



## P-2150

### OVERALL SIZE

- 30"W x 19 $\frac{3}{4}$ "D x 75"H

### COMPARTMENT SIZE

- 3.7 cu. ft. fluid compartment
- 24 $\frac{1}{4}$ "W x 18"D x 15"H
- 9 cu. ft. blanket compartment
- 24 $\frac{1}{4}$ "W x 18"D x 36"H

- **tri▼temp Range:**
  - Blankets: 90°F - 160°F
  - Irrigation: 90°F - 150°F
  - Injection: 90°F - 104°F
- **Total Capacity:**
  - 24 1-liter bags
  - 20 1-liter bottles
  - 9-12 blankets (upper chamber)
  - 24-30 blankets (lower chamber)
- Separate **tri▼temp** temperature controls for each compartment
- Pass Through Option Available
- Data Logger



# warming cabinets

## blanket

### features

Elite Series blanket warmers feature windowed door, programmable control with temperature lock and four digit LED display and single pane energy efficient coated glass window door allows inventory observation.

### P-2055

#### OVERALL SIZE

32"W x 29¼"D x 71"H

#### COMPARTMENT SIZE

- 10.3 cubic ft. upper
- 26½"W x 23¾"D x 28⅝"H
- 10.3 cubic ft. lower
- 26½"W x 23¾"D x 28⅝"H

- Temp Range:
  - Blankets: 98°F - 160°F
- Total Capacity: 60-70 blankets
- Auto time programming
- One compartment, two separate doors
- No fans or lint traps
- 5" casters
- Pass through option available



Heated center shelf



### P-2210

#### OVERALL SIZE

- 18½"W x 22"D x 22"H

#### COMPARTMENT SIZE

- 2.5 cubic ft.
- 15½"W x 18¾"D x 15¼"H

- Temp Range:
  - Blankets: 90°F - 160°F
- A warming shut-off system, separate from the electronic control, prevents overheating
- Interior blue lighting



### P-2220

#### OVERALL SIZE

- 18½"W x 22"D x 28"H

#### COMPARTMENT SIZE

- 3.5 cubic ft.
- 15½"W x 18¾"D x 21¼"H\*

\* (9 9/16"H above & below the center shelf)

- Temp Range:
  - Blankets: 90°F - 160°F
- A warming shut-off system, separate from the electronic control, prevents overheating
- Interior blue lighting



### P-2230

#### OVERALL SIZE

- 24"W x 24½"D x 22"H

#### COMPARTMENT SIZE

- 4 cubic ft.
- 21"W x 21½"D x 15¼"H

- Temp Range:
  - Blankets: 90°F - 160°F
- A warming shut-off system, separate from the electronic control, prevents overheating
- Interior blue lighting

Heated center shelf



### P-2240

#### OVERALL SIZE

- 24"W x 24¼"D x 35"H

#### COMPARTMENT SIZE

- 7.5 cubic ft.
- 21"W x 21¾"D x 28¼"H\*

\* (13 1/8"H above & below the center shelf)

- Temp Range:
  - Blankets: 90°F - 160°F
- A warming shut-off system, separate from the electronic control, prevents overheating
- Interior blue lighting



### P-2112

#### OVERALL SIZE

18½"W x 23½"D x 22"H

#### COMPARTMENT SIZE

- 2.5 cubic ft.
- 15½"W x 18¾"D x 15¼"H

- Temp Range:
  - Irrigation: 90°F - 150°F
  - Injection: 90°F - 104°F
- Total Capacity:
  - 16 1-liter bags or
  - 12 1-liter bottles
- Glide out basket
- Casters (option)



### P-2122

#### OVERALL SIZE

18½"W x 23½"D x 28"H

#### COMPARTMENT SIZE

- 3.5 cubic ft.
- 15½"W x 18¾"D x 21¼"H

- Temp Range:
  - Irrigation: 90°F - 150°F
  - Injection: 90°F - 104°F
- Total Capacity:
  - 24 1-liter bags or
  - 24 1-liter bottles
- Glide out basket
- Casters (option)

## warming cabinets

fluid

### features

Fluid warmers feature 4 digit LED display, window doors, programmable temperature control and interior lighting.



### P-2030-TC

#### OVERALL SIZE

26"W x 33"D x 28"H

### features

- Stainless steel
- 4" stainless steel casters with brakes
- 400 lbs. total weight capacity at 200 lbs. per shelf
- Top shelf with 3-sided guard rail
- Bottom shelf with 4-sided rail.
- Compatible with the P-2032, P-2230 and P-2240 warming cabinets

shown with the  
P-2010-S warmer



### P-2010-TC

#### OVERALL SIZE

23¼"W x 32¼"D x 43"H

### features

- Stainless steel
- 4" casters with brake & swivel lock
- 300 lbs. total weight capacity at 100 lbs. per shelf
- Shelf Dimensions  
17"W x 27"D x 12¼"H
- Compatible with the P-2010-S warming cabinet



### P-2010-CB

#### OVERALL SIZE

21¼"W x 29¾"D x 6<sup>13</sup>/<sub>16</sub>"H

### features

- Platform
- 4" stainless steel casters with brakes
- 350 lbs. total weight capacity
- Perimeter bumper
- Cutouts in platform capture warming cabinet feet.
- Compatible with the P-2010-S warming cabinet

# warming cabinet features



## Fluid Warmers

- **Radiant Heat Technology**
- Safety coated interior
- Fully insulated
- 120 VAC
- Heavy duty hardware
- Glass door
- Perimeter bumper (P-2130 only)

### Safety Features

- Built-in speaker for audible feedback
- Safety override feature prevents overheating
- Lockout feature

### Options

- Optional left hand door swing available
- Option of casters or fixed legs

### Warranty

- One Year Service
- Three Years on Parts
- Lifetime Technical Support



### Models

- P-2112
  - P-2122
  - P-2130
  - P-2140\*
- \* combination fluid/blanket warmer



## Blanket Warmers

- **Radiant Heat Technology**
- Safety coated interior
- Fully insulated
- 120 VAC
- Heavy duty hardware
- Glass door (P-2010-S solid door)

### Safety Features

- Built-in speaker for audible feedback
- Safety override feature prevents overheating
- Lockout feature

### Options

- Optional left hand door swing available
- Option of casters or fixed legs
- Perimeter bumper is optional (excludes P-2055)

### Warranty

- One Year Service
- Three Years on Parts
- Lifetime Technical Support



### Models

- P-2010-S
- P-2012
- P-2022
- P-2032
- P-2040
- P-2055
- P-2210
- P-2220
- P-2230
- P-2240



## Blanket & Fluid Warmers

- **Radiant Heat Technology**
- **tri▼temp**
- Safety coated interior
- Fully insulated
- 120 VAC
- Heavy duty hardware
- Glass or solid door

### Safety Features

- Built-in speaker for audible feedback
- Safety override feature prevents overheating (Deluxe models only)
- Lockout feature (Deluxe models only)

### Options

- Optional left hand door swing available
- Option of casters or fixed legs (excludes P-2150)
- Optional pass through warmer (P-2150 only).

### Warranty

- One Year Service
- Three Years on Parts
- Lifetime Technical Support



### Models

- P-2147
- P-2150
- P-2420
- P-2430





# ivNow

## individual warming pods

### Modular Configuration

Wall mount:

One to six (1-3) bag capacity

Countertop:

One to three (1-3) bag capacity

Pole mount:

One to three (1-3) bag capacity

### Temperature:

- Heats to 104°F (40°C) in 30 minutes or less (slightly longer heating time is required for 3L bags)
- Accuracy of +0°/-2°C (+0°/-3°F)

### Electronic Control

Monitors the length of time the bag has been held at temperature, displayed additionally by a status button. The control alerts users when a fluid bag has been held at temperature longer than 14 days.

### Warranty

1 year parts and labor



### features

The ivNow® fluid warmer quickly warms and maintains the temperature of injection/intravenous and irrigation solutions prior to their use. The specially contoured warming module cradles solution bags in 0.5-, 1-, 2- & 3-liter sizes. ivNow® units can be customized in the field by stacking or removing pods from the configuration. Units can be used for countertop use or wall mounted with the provided mounting plate.



### ivNow-1

The ivNow-1 may be placed directly on the countertop or fitted with an optional bracket for mounting on a stand or wall.

P-2301-A

#### Mount

- Equipment pole stand
- Wall mounting brackets

Capacity

One (1) bag, up to 3l \*



### ivNow-2

The ivNow-2 may be placed directly on the countertop or fitted with an bracket for mounting on a stand or wall.

P-2302-A

#### Mount

- Equipment pole stand
- Wall mounting brackets

Capacity

Two (2) bags, up to 3l \*



### ivNow-3

The ivNow-3 may be placed directly on the countertop or fitted with an bracket for mounting on a stand or wall.

P-2303-A

#### Mount

- Equipment pole stand
- Wall mounting brackets

Capacity

Three (3) bags, up to 3l \*