


Galaxy Gen 3 4-Pod Mobile Stand Light



Touch Control Pad Detail

Level Indicator
Dimming Control
5 Levels
Power
LED Medical Lighting
StarTrol.com

<u>Galaxy Gen3 Mobile 4 Pod Specifications</u>	
Model Number	4x4-GSM
Input Voltage	100-240 VAC
Input Current	0.33 Amps @ 120 VAC
Input Wattage	40 Watts
Light Output	70,000 Lux @ 1 Meter
Dimming	5 Levels
Weight	53 Lbs.
Power Cord Length	10'
Max Head Height	80"
Min Head Height	41"