

Generation 3 Galaxy 4X4

Mobile Exam Light

Part #: 4X4-GSM

Introduction

The **StarTrol** Galaxy Generation 3, the newest member of our lighting family. The LED's are brighter than our current lights and features five level dimming. The new sleek design, features a fully enclosed light head for improved clean ability. *StarTrol* is proudly made in the United States, with an unmatched 10 year limited warranty.

Standard Features

Touch Sense on / off with 5 step dimming, along with level indicator lights Multiple configurations: Mobile, Wall, Single, and Dual Ceilings feature spring arm suspension Removable anodized aluminum handle Enclosed light head for improved clean ability

10 year limited warranty

Architectural Specifications StarTrol Gen 3 4X4-GSM has an 8.75" diameter aluminum light head featuring 16 LEDs (70,000 lux at a meter). This light has an arm reach of 31". The light head has a 700° rotation at the arm and head pivot point and 330° at the light head joint. The overall height of the light is 60.5", the dimensions of mobile base are 21" x 20.1" x 12.5". Electrical power is supplied by a 10' hospital grade power cord (included).

Technical Data

Part number: 4X4-GSM

Input voltage: 100-240 VAC 50-60Hz Auto Adjusting

Input Current: .33 Amps @ 120 VAC Light Output: 70,000 lux at a meter

Spot Size: 8 inchés

Color Temperature: 4500K Color rendering index: 95

R9: 99

LED Life: 50,000+ hours **Number of LEDs: 16**

Dimming: 5 Level Touch Sense

Total Reach: 31 inches

Base dimension: 21" X 20.1" 12.5"



Accessories

Additional Aluminum handle (Part# HDN) Replacement Power cord (Part# PWC)