




# Typical Drop Ceiling Install Galaxy Gen 3 Single 4-Pod Light



## Touch Control Pad Detail 5 Level Dimming



### Galaxy Gen3 4-Pod Single Specifications

Model Number	4x4-GS1
Input Voltage	100-240 VAC
Input Current	0.33 Amps @ 120 VAC
Input Wattage	40 Watts
Light Output	70,000 Lux @ 1 Meter
Moment Load	70 ft-lbs
Max Weight	49 Lbs.
Min. Ceiling Height	93" @ 80" Pivot Hgt.
Max. Ceiling Height	148" @ 80" Pivot Hgt.

### Notes:

1. Customer is responsible for consulting with structural engineer to determine design of support structure.
2. Electrical connections to be completed by licensed electrical contractor in compliance with all local codes.
3. Ceiling Box can also be mounted at finished ceiling height with additional ceiling cover.

## Mounting Info

