



Generation 3 Galaxy 8X4

Mobile Light
Part #: 8X4-GSM

Introduction

The **StarTrol** Galaxy Generation 3, the newest member of our lighting family. The LED's are brighter than our current lights and features five level dimming. The new sleek design, features a fully enclosed light head for improved clean ability. **StarTrol** is proudly made in the United States, with an unmatched 10 year limited warranty.

Standard Features

Touch Sense on / off with 5 step dimming, along with level indicator lights
Multiple configurations: Mobile, Wall, Single, and Dual
Ceilings feature spring arm suspension
Removable anodized aluminum handle
Enclosed light head for improved clean ability
10 year limited warranty

Architectural Specifications

StarTrol Gen 3 8X4-GSM has a 12.5" diameter aluminum light head featuring 32 LEDs (140,000 lux at a meter). This light has an arm reach of 31". The light head has a 700° rotation at the arm and head pivot point and 330° at the light head joint. The overall height of the light is 60.5", the dimensions of mobile base are 21" x 20.1" x 12.5". Electrical power is supplied by a 10' hospital grade power cord (included).

Technical Data

Part number: 8X4-GSM
Input voltage: 100-240 VAC 50-60Hz Auto Adjusting
Input Current: 0.67 Amps @ 120 VAC
Light Output: 140,000 lux at a meter
Spot Size: 10 inches
Color Temperature: 4500K
Color rendering index: 95
R9: 99
LED Life: 50,000+ hours
Number of LEDs: 32
Dimming: 5 Level Touch Sense
Total Reach: 31 inches
Base dimension: 21" x 20.1" x 12.5"



Accessories

Additional Aluminum handle (Part# HDN)
Replacement Power cord (Part# PWC)
Lens (Part# LNS)