



# Galaxy Gen 3 Wall Mounted 8-Pod Light



## Touch Control Pad Detail 5 Level Dimming



### Notes:

1. Customer is responsible for consulting with structural engineer to determine structural backing and hardware.
2. Electrical connections to be completed by licensed electrical contractor in compliance with all local codes.
3. Light can be equipped with cord or hard wired.

## Galaxy Gen3 8- Pod Specifications

Model Number	8x4-GSW
Input Voltage	100-240 VAC
Input Current	0.67 Amps @ 120 VAC
Input Wattage	80 Watts
Light Output	140,000 Lux @ 1 Meter
Dimming	5 Level Touch Sense
Weight	49 Lbs.

## Mounting Bracket

Typical Mounting Dimensions for 80" Pivot Height

