



Generation 3 Galaxy 8X4

Wall Mounted Light
Part #: 8X4-GSW



Introduction

The **StarTrol** Galaxy Generation 3, the newest member of our lighting family. The LED's are brighter than our current lights and features five level dimming. The new sleek design, features a fully enclosed light head for improved clean ability. **StarTrol** is proudly made in the United States, with an unmatched 10 year limited warranty.

Standard Features

- Unlimited Rotation at the main pivot point
- Robust spring arm- that never needs adjustment
- Touch Sense on / off with 5 step dimming, along with level indicator lights
- Multiple configurations (Wall, Single, and Dual Ceilings)
- Removable anodized aluminum handle
- Enclosed light head for improved clean ability
- 10 year limited warranty

Accessories

- Additional aluminum handle (Part #: HDN)
- Replacement power cord (Part #: PWC)
- Lens (Part #: LNS)

Architectural Specifications

StarTrol Gen 3 8X4 GSW has an 12.5" diameter aluminum light head featuring 32 LEDs (140,000 lux at a meter). This light has an arm reach of 52" to 76" in 3" increments. with unlimited rotation at the main pivot point. The light head has a 700° rotation at the arm and head pivot point and 330° at the light head joint. Electrical power is supplied via a 10' hospital grade power cord (included), can also be hardwired.

Technical Data

- Part number:** 8X4-GSW
- Input voltage:** 100-240 VAC 50-60Hz Auto Adjusting
- Input Current:** 0.67 Amps @ 120 VAC
- Light Output:** 140,000 lux at a meter
- Spot Size:** 10 inches
- Color Temperature:** 4500K
- Color rendering index:** 95
- R9:** 99
- LED Life:** 50,000+ hours
- Number of LEDs:** 32
- Dimming:** 5 Level Touch Sense
- Total Reach:** 60" standard reach, Customizable from 52" to 76" in 3" increments.