



basin stands

and kick buckets



features

- Elevated ring basin bracket
- Mobile and easy to clean
- Open frame design allows access from all sides
- 2" diameter double ball bearing swivel casters



P-78

single basin stand

Near chrome powder coat finish
13³/₄"D x 14"W x 33³/₄"H

- One removable 13¹/₂" diameter, 7 quart stainless steel basin
- Steel construction

P-1078-SS

single basin stand

Stainless steel
13³/₄"D x 14"W x 33³/₄"H



P-1078-W/S-SS

single basin stand with shelf

Stainless steel
13³/₄"D x 14"W x 33³/₄"H

- One removable 13¹/₂" diameter, 7 quart stainless steel basin
- Stainless steel construction

P-1079-W/S-SS

double basin stand with shelf

Stainless steel
14"D x 28"W x 33³/₄"H



P-79

double basin stand

Near chrome powder coat finish
13³/₄"D x 28"W x 33³/₄"H

- Two removable 13¹/₂" diameter, 7 quart stainless steel basins
- Steel construction

P-1079-SS

double basin stand

Stainless steel
14"D x 28"W x 33³/₄"H

features

- Balanced, tip-resistant
- Perimeter bumper to protect walls
- Electro-polished, stainless steel frame
- 2" diameter double ball bearing swivel casters for easy mobility

P-1020-SS

kick bucket

Stainless steel
Overall size
13³/₄"diameter x 13¹/₄"H
• Removable 12 quart, stainless steel bucket



P-1025-SS

sponge receptacle

Stainless steel
Overall size
16¹/₂"diameter x 7¹/₂"H
• Removable 8¹/₂ quart, stainless steel basin

