


Datascope Trio™



Perhaps the most clever compact yet.

Datascope Trio™

We've raised the standard on what's standard.

Designed specifically to address the needs of today's highly diverse patient population,

Trio packs all the features you need in a compact, cost-effective solution. But it doesn't end there. Trio raises the standard in compact patient monitoring by looking beyond bare necessity, adding features that make a good monitor a great one. It successfully partners a high-resolution, easy-to-view color display with a truly practical, lightweight design. Trio also offers valuable options to further extend your capabilities. The result: a highly adaptable, portable package that you'll want to use wherever possible – from bedside to transport.

[BASIC SYSTEM]

- 8.4" TFT color display with 800 x 600 resolution
- 4-trace display with ECG cascade
- 3 or 5-lead ECG
- ECG-derived Impedance Respiration
- Non-Invasive Blood Pressure
- Masimo SET® motion tolerant SpO₂
- Continuous Temperature
- Full Graphic and List Trends of all monitored parameters
- Front panel keypad featuring Navigator Control Knob and Quick Action Keys
- Conveniently located Alarm Light
- Connectivity to other devices
- Built-in power supply for direct AC connection
- Patented design with all-in-one fold-away handle and bed rail hook
- Customizable default settings for each patient size
- Lightweight, compact design

[AVAILABLE OPTIONS]

- Nellcor® OxiMax® SpO₂
- Invasive Blood Pressure
- Lithium Ion Battery
- Integral, Dual-Trace Recorder


A HISTORY OF INNOVATION



1964
Carditron
First cardiac monitor with synchronized ECG



1986
Accusat
Datascope's first stand-alone pulse oximeter



1989
Accutorr™
All-in-one blood pressure monitor with rotating screen



1992
Passport™ LT
Reduced parameter version of the popular Passport monitor



1996
Accutorr™ Plus
Accutorr reliability, plus advanced feature-set and compact design



2001
Radical™
Stand-alone pulse oximeter offering Masimo SET technology



2003
Trio™
Truly compact, portable bedside monitor

Added flexibility without the added cost.

Through exceptional versatility, Trio™ delivers exceptional value.

With its innovative design, Trio meets the demands of today's transient and highly mobile patients. Whether supporting the short-term needs of an outpatient, the immediate needs of a transport patient, or longer term needs at the bedside, Trio is up to the task.

Mounted to a wall or hooked to a bed rail, Trio's well-organized display makes data analysis a pleasure. The brilliant, 8.4" color display provides exceptional viewing, no matter what the setting. Trio's unique, all-in-one fold-away handle and bed rail hook adds a new dimension to compact monitoring – ultimate flexibility. Just unplug the power cord and Trio automatically switches to battery power for uninterrupted, in-transit monitoring.

All this, packed into a durable, lightweight monitor made to pick up and go.



Trio makes transport easy

with a built-in, fold-away handle and bed rail hook that mounts directly to the bed. Add a rolling stand or wall mount with adjustable arm and go from transport to bedside instantly.





A lightweight and compact design
provides ease of portability, allowing
you to bring Trio effortlessly from one
place to the next.

You may call it clever, innovative,
remarkably versatile. We just call it Trio™.

With Trio, select only those options critical to you.

A dual-trace, integral recorder with one-touch printing allows you to perform a variety of functions such as print on alarm, trend print, and continuous real-time print. Trio satisfies your requirement for extended battery life with its Lithium Ion technology. Trio also offers a choice of invasive blood pressure, Masimo SET® motion tolerant SpO₂ or Nellcor® OxiMax® SpO₂.

Add to all this a host of mounting alternatives, and you have a versatile, affordable compact monitor uniquely adaptable to your demanding environment.



Masimo SET

Masimo SET®
Motion Tolerant
SpO₂



Dual-Trace
Recorder



Lithium Ion Battery



- 1 **ETHERNET CONNECTOR** for quick and simple software upgrades.
- 2 **ANALOG OUTPUT** for ECG analog signal output to other devices.
- 3 **SERIAL OUTPUT** for connectivity to other devices.
- 4 **BUILT-IN POWER SUPPLY** for direct AC connection.



Duo™



Accutorr® V



Trio™



Passport® V



Spectrum®



Spectrum® OR



AS 3000™

A monitor you can trust, from a company you can trust.

Our service organization is one of the largest service teams in the industry. With our selection of highly personalized service contracts, parts and accessories, loaner and rental programs, we can protect and add value to your equipment investment. Multi-level training is available for biomedical engineers at your site or ours. We offer standard and customized programs designed to help you reduce service costs and maximize up-time.

Visit mindray.com to learn more about our full line of monitoring solutions, including our wide variety of parts and accessories.

Mindray DS USA, Inc.
800 MacArthur Boulevard
Mahwah, NJ 07430 USA
Tel 1 800 288 2121
or 201 995 8000
Fax 1 800 926 4275

www.mindray.com

mindray™
NORTH AMERICA

Duo™, Trio™ and AS3000™ are trademarks of Mindray DS USA, Inc. Accutorr® is a trademark of Mindray DS USA, Inc. and is registered in the U.S. and other countries. Spectrum® and Passport® are U.S. registered trademarks or trademarks of Mindray DS USA, Inc. Masimo SET® is a U.S. registered trademark and Radical™ is a U.S. trademark of Masimo Corp. Nellcor® and OxilMax® are U.S. registered trademarks of Nellcor Puritan Bennett Inc. Datascope is a trademark or registered trademark of Datascope Corp. or its subsidiaries and is used under license by Mindray DS USA, Inc.



© 2010 Mindray DS USA, Inc.
0002-08-13003 Rev E