

Maximum Power.

Maximum Performance.

Proven Quality.

- Unique, straight-forward user interface – a hallmark of the Spectrum® series – easily configurable to meet specific requirements
- Bright, easy-to-read, user selectable real-time O₂, CO₂, N₂O, and Agent waveform and numerics display, status messages, and alarms greatly simplify procedural monitoring – even at a distance
- Breath-by-breath O₂, CO₂, and N₂O monitoring with automatic 5-agent identification. Provides for real-time display of all critical inspired and expired gases
- Compact, light-weight design ensures convenient mounting on anesthesia delivery systems
- Spectrum OR provides comprehensive vital signs monitoring capability with optional remote color display available for enhanced visibility

Gas Module 3™

State-of-the-art comprehensive vital signs monitoring *and* multi-gas analysis for any surgical environment.



**Spectrum OR and
Gas Module 3**

We proudly offer the Gas Module 3, the most powerful performance extension to our Spectrum OR line of monitors – already a standard of vital signs monitoring excellence worldwide. The Gas Module 3 delivers state-of-the-art gas monitoring and analysis capabilities: identification and measurement of inspired and expired O₂, CO₂, N₂O, and 5 anesthetic agents – Isoflurane, Sevoflurane, Desflurane, Halothane, and Enflurane. Providing impressive capabilities and highly reliable performance, this leading-edge breath-by-breath multi-gas analyzer was designed to meet the comprehensive anesthesia monitoring requirements of virtually every hospital – whatever its size, specialty, or patient base.

Gas Module 3™

So accurate. So reliable. So simple to use.

Our compact, lightweight monitors are perfectly suited to the surgical setting. The simplified cable management system ensures fast, convenient placement on the anesthesia delivery system.



Large, bright displays, easily visible from a distance, and a simple user interface make this vital signs/multi-gas analysis solution the logical choice.

Mindray DS USA, Inc.
800 MacArthur Boulevard
Mahwah, NJ 07430
Tel 1 800 288 2121
Tel 1 201 995 8000
Fax 1 800 926 4275
www.mindray.com

Gas Module 3™ is a trademark of Mindray DS USA, Inc.
Spectrum™ is a U.S. registered trademark or trademark of Mindray DS USA, Inc.
Mindray™ is a trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Datascope is a trademark or registered trademark of Datascope Corp. or its subsidiaries
and is used under license by Mindray DS USA, Inc.
Subject to change. © 2010 Mindray DS USA, Inc.

P/N: 0002-08-6010 Rev C