


MEDI CAPTURE

Medical Imaging Made Easy™



MediCap® USB200

The MediCap USB200 is the medical industry's easiest full-featured video recording solution. Your videos and still images are recorded directly to a USB flash drive or optional hard drive. The recorded videos are DVD-quality, but without all the hassles of DVD disks.



Medical Video Recorder

Records from your endoscope, arthroscope, surgical camera, ultrasound, etc. Saves to a flash drive or hard drive that can be accessed from any computer on your network.



Intuitive Use

It's so easy, no training is needed. Just plug it in and go. You'll be recording videos within minutes of opening the box.



Medical Grade

With its steel case and fluid resistant front panel, the USB200 is built to withstand the rigors of the OR, and it's so compact it will fit on any cart.

Powerful Options

compatible – versatile – affordable



USB Hard Drive

Perfect for archiving
100+ hours of video
on a single device



USB Keyboard

Input patient data to
organize your photos
and videos



Network Adapter

Turn your USB200
into a mapped drive
on your network



Foot Pedal

Connect your USB200
to your scope's trigger
or to a foot pedal



Also Available...



MediCap® USB300

Step up to high definition. The USB300 has all the same great features as the USB200, plus HD recording — at a surprisingly affordable price.



MediCap® USB200 Specs

| | |
|-------------------------------|---|
| Device Type | Medical Video Recorder Records video and still images from standard medical video sources to a USB flash drive or optional USB hard drive |
| Case | Rugged metal case with high-impact plastic faceplate |
| LCD Monitor | 2 1/2", color (for video confirmation only, not for diagnosis) |
| Speaker | 200Hz-15kHz, 1W, (for audio confirmation only, not for diagnosis) |
| Buttons | Sealed membrane, fluid resistant |
| Video In/Out | PAL or NTSC format using S-video or composite BNC |
| Audio In/Out | Standard 3.5mm stereo phone jacks |
| USB Support | USB2.0, high speed, 1 front connector, 3 rear connectors |
| Networking | Ethernet T10/100 connection via Network Upgrade Kit (optional) |
| Footswitch/Trigger | Hands-free capture Standard 3.5mm mini jack |
| Mechanical | Size: 240 x 210 x 63mm / 9.5 x 8.27 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.) |
| Environmental | Storage: -40° to +85° C Operating: -20° to +40° C |
| Power | AC 100-240V 50/60Hz 3-prong AC jack |
| Media | USB flash drives or external USB hard drive |
| Image/Video Resolution | NTSC: 720 x 480 pixels PAL: 720 x 576 pixels |
| Image Formats | JPEG, TIFF, BMP, DICOM |
| Video Format | H.264 (DVD-quality) |
| Audio Format | AAC, 2 channels |
| Certifications | CE 93/42/EEC, IEC 60601-1:2005, ANSI/AAMI ES60601-1:2005, CAN/CSA C22.2 No.60601-1:2008, E361228, FCC Part 15, FDA Class 1, RoHS, PSE, GOST-R |
| Kit Contents | MediCap USB200, USB flash drive, video cables, power cable, user guide |

MediCapture® Inc.

2250 Hickory Road, Suite 200
Plymouth Meeting, PA 19462 USA

(888) 922-7887
(503) 445-6935

www.medicapture.com
sales@medicapture.com

